**20 GENESYS<sup>9</sup> SPECTROPHOTOMETER**  
**Operator's Manual**  
**Manual del Operador**  
**Manuel de l'opérateur**  
**Benutzerhandbuch**

**Thermo** Spectronic

Copyright © 2000, Spectronic Unicam  
All rights reserved.

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
**General safety notes  
used in this manual**

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
This symbol alerts you to important information about using the instrument. Be sure to read and follow the associated instructions carefully.

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
This symbol alerts you to potential electrical hazards. Be sure that only qualified persons perform the related procedures.

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
This symbol alerts you to hot surfaces. Be sure to read and follow the associated instructions carefully.

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This symbol alerts you to the presence of static-sensitive devices. Be sure to read and follow the associated instructions carefully.

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**LIMITED WARRANTY**

Spectronic Unicam warrants the 20 GENESYS9 against

defects in material and workmanship for a period of (3) years from the date of delivery. Related accessories are warranted against defects in material and workmanship for a period of (1) year. For warranty repair the 20 GENESYS9 or accessory will be returned to a service depot which is authorized by Spectronic Unicam to perform warranty repair and service. The prompt return of your warranty registration card can assist Spectronic Unicam in expediting warranty repair should a problem occur.

This warranty covers all parts (except those specified below) and applies only to equipment which has been installed and operated in accordance with the operator's instruction manual and which has been serviced only by authorized Spectronic Unicam=s distributors or service personnel. This warranty does not apply to equipment and accessories that have been modified or tampered with in any way, misused, or damaged by accident, neglect, or conditions beyond Spectronic Unicam=s control.

This warranty does not apply to lamps, glassware, and similar expendable components. However, such parts and components may be warranted by their manufacturer.

Spectronic Unicam is not responsible under this warranty for loss in operating performance due to environmental conditions.

**THIS WARRANTY IS IN LIEU OF ALL WARRANTIES EXPRESS, IMPLIED, OR STATUTORY, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF FITNESS FOR A PARTICULAR PURPOSE OR MERCHANTABILITY OR OTHERWISE,** and states Spectronic Unicam=s entire and exclusive liability and the customer's exclusive remedy for any claim in connection with the sale or furnishing of services, goods, or parts, their design, suitability for use, installation, or operations. SpectronicUnicam will in no event be liable for any direct, indirect, special, or consequential damages whatsoever, including loss of goodwill, whether grounded in tort (including negligence), strict liability or contract, and Spectronic Unicam=s liability under no circumstances will exceed the contract price for the goods and/or services for which liability is claimed.

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## Table of Contents

|   |    |
|---|----|
| Description.....                                    | 1  |
| Components.....                                     | 2  |
| Keyboard layout.....                                | 4  |
| Specifications.....                                 | 6  |
| Environmental and electrical requirements.....      | 8  |
| Basic operation.....                                | 10 |
| Setting up the instrument.....                      | 12 |
| Connecting to external devices.....                 | 14 |
| Starting up the instrument.....                     | 14 |
| Absorbance and %transmittance<br>measurements.....  | 16 |
| Concentration measurements<br>using a factor.....   | 18 |
| Concentration measurements<br>using a standard..... | 20 |
| Setting absorbance and %transmittance.....          | 22 |
| Printing your results.....                          | 22 |
| Using the utility functions.....                    | 24 |
| Concentration units.....                            | 24 |
| Language.....                                       | 26 |
| RS232.....  | 26 |
| Paper feed.....                                     | 26 |
| Printer setup.....                                  | 28 |
| Lamp alignment.....                                 | 28 |
| Hours of lamp operation.....                        | 30 |
| Initial wavelength.....                             | 32 |
| Beeper control.....                                 | 32 |
| Standby.....  | 34 |
| Software revision.....                              | 34 |
| Accessories.....                                    | 36 |
| Removing the cell holder.....                       | 38 |
| Installing the cell holder.....                     | 38 |
| Maintenance & troubleshooting.....                  | 40 |
| Routine care.....                                   | 40 |
| Changing the lamp.....                              | 42 |
| Changing the fuse.....                              | 48 |
| Replacement parts.....                              | 52 |
| Messages.....                                       | 52 |

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#### **NOTE**

This operator's manual contains information, instructions and specifications for the 20 GENESYS<sup>9</sup> spectrophotometer that were believed accurate at the time this manual was written. However, as part of Spectronic Unicam's on-going program of product development, the specifications and operating instructions may be changed from time to time. Spectronic Unicam reserves the right to change such operating instructions and specifications. Under no circumstances shall Spectronic Unicam be obligated to notify purchasers of any future changes in either this or any other instructions or specifications relating to Spectronic Unicam's products, nor shall Spectronic Unicam be liable in any way for its failure to notify purchasers of such changes.

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### Description

The 20 GENESYS9 spectrophotometer is an easy-to-use instrument that performs absorbance, %transmittance and concentration measurements within the wavelength range of 325 to 1100 nanometers. Its features include:

- X Two-line, 20-character LCD display
- X Tactile keyboard that clicks to indicate a key was pressed
- X Optional cell holders available for a wide range of cuvettes, test tubes and longpath cells
- X Optional internal printer
- X Built-in RS232 interface

**Note:** *If you plan to use the RS232 interface, please request a copy of the RS232 command set. In the United States, you can call 1-800-654-9955; in other countries, contact your distributor.*

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**Components**

1. On/Off switch
2. LCD display
3. Sample compartment door
4. Keyboard
5. Optional built-in printer
6. Lamp compartment door

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### Keyboard layout

1. Display - 20-character, 2-line LCD
2. Soft key 1 - Function varies depending on screen; generally Escape, Back up or Clear
3. Soft key 2 - Function varies depending on screen; generally Enter, Accept or Continue
4. Scroll keys - Used to scroll through menus and enter numeric values
5. Wavelength controls - Increase and decrease the wavelength setting
6. 0 Abs/100%T - Automatically sets the instrument to zero absorbance (100%T)
7. A/T/C - Switches between absorbance, %transmittance and concentration modes
8. Utility - Accesses instrument set-up, diagnostics and other functions
9. Print - Sends currently displayed data to selected printer

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### Specifications

|                          |  |
|--------------------------|--|
| Wavelength range         | 325 to 1100 nm   |
| Bandpass                 | #8 mm  |
| Stability                | #3 mA/hour drift   |
| Stray light              | #0.1%T, when measured at 340 and 400 nm                                    |
| Wavelength repeatability | ∓0.5 nm  |
| Wavelength accuracy      | ∓2.0 nm  |
| Noise (@ 500 nm)         | #1 mA at 0A and 2 mA at 2A, peak-to-peak (15 seconds)                      |
| Photometric accuracy     |  |
| 0.0 to 0.3 A             | 0.003 A  |
| 0.301 A to 2.5 A         | 1.0%   |
| Photometric range        | 0-125%T, -0.1-2.5A, 0-1999C  |
| Light source lifetime    | Visible: ~1000 hours   |
| Grating                  | 1200 lines/mm  |
| Data output              | ∃ Two-line, 20-character LCD display<br>∃ RS232C port<br>∃ Centronics port |
| Dimensions               | 30 cm (12") W x 33 cm (13") D x 19 cm (7") H                               |
| Power requirements       | Selected automatically; 100 to 240 Volts; 50 to 60 Hz, 1.0 Amp             |

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### Environmental and electrical requirements

Your 20 GENESYS9 spectrophotometer has been designed to operate under the environmental and electrical requirements listed below.

#### Line voltages

100 - 240V  $\nabla$  10% (max.)  
50 - 60 Hz  $\nabla$  10%

**Operating environment** - The instrument meets the specifications on the previous pages under the following conditions after a 30-minute warm-up period.

| Ambient temperature          | Relative humidity |
|------------------------------|-------------------|
| 50EF to 75EF (15EC to 24EC)  | 20% to 80%        |
| 76EF to 85EF (25EC to 29EC)  | 20% to 70%        |
| 86EF to 95EF (30EC to 35EC)  | 20% to 60%        |
| 96EF to 105EF (36EC to 40EC) | 20% to 50%        |

Temperature should be maintained at  $\nabla$ 4EF ( $\nabla$ 2EC).  
Relative humidity should be maintained to  $\nabla$ 5%.

#### Storage environment

-40EF to 140EF (-40EC to 60EC) - Relative humidity not to exceed 60%. Allow instrument to adjust to room temperature for 24 hours after taking it out of storage.

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**Installationskategorie II**

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### Basic operation

Using the 20 GENESYS9 spectrometer, you can perform absorbance and %transmittance measurements, and determine concentrations using either a standard solution or a conversion factor. Regardless of the type of measurement you want to perform, you follow similar steps:

- X Select the mode (A, %T, C) and wavelength
- X Measure the blank
- X Enter the standard value or factor (concentration mode only)
- X Measure the samples

You may also change default settings for the instrument and printer by accessing the SETUP function. Setup includes language, RS232 parameters, printer format, lamp alignment and usage, beeper control, standby mode and others.

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## Setting up the instrument

Unpack the shipping carton and verify that you have received all the items listed on the checklist. If any item is missing, contact your local distributor.

Place the instrument on a flat, even surface that is:

- X As far as possible from any strong electric or magnetic fields and from any electrical device that may generate high-frequency fields
  - X Free of dust, corrosive gases and strong vibrations
- Remove any obstructions or materials that could hinder the flow of air under, behind and around the instrument. Connect the female end of the power cord into the A/C power connector on the back of the instrument. Plug the power cord into a grounded outlet with the appropriate voltage. Turn the power switch to ON (\*=ON, ±=OFF).

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### Connecting to external devices

Your 20 GENESYS9 spectrometer supports a variety of external devices:

- X RS-232 devices - Supports any RS-232 device such as a terminal, computer or printer. The instrument behaves as a DCE device.
- X Printers - Supports Centronics9-compatible devices

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### Starting up the instrument

When you turn on your 20 GENESYS9 spectrometer, it performs its power-on sequence. This sequence includes checking the software revision, initializing the filter wheel and the monochromator. The power-up sequence takes about two minutes to complete. Allow the instrument to warm up for 30 minutes before using it.

**Note:** *Be sure that the cell holder is empty before turning on the instrument.*

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## Concentration measurements using a standard

**1**

Press **A/T/C** to select the concentration mode. The current mode appears on the display.

**2**


Press **nm**  $\square$  or **nm**  $\square$  to select the wavelength.

**Note:**  *Holding either key will cause the wavelength to change more quickly.*

**3**

Insert your blank into the cell holder and close the sample door.

**Note:**  *Position the cell so the light (indicated by arrow in drawing) passes through the clear walls.*

**4**

Press **0 ABS/100%T** to zero the blank, then remove the blank and insert your standard.


**5**

Press the **Standard** soft key and use the  $\square$  and  $\square$  keys to enter the concentration of the standard, then press the **Set C** soft key to calculate and display the factor for the selected standard.

**6**

Remove the standard and insert your sample into the cell holder. The calculated sample concentration appears on the LCD display.

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

### Setting absorbance and %transmittance

In addition to taking absorbance and %transmittance measurements, you can also set specific values for the absorbance or %transmittance on the display.

- X Press the  $\square$  and  $\square$  keys to set the absorbance or %transmittance value on the display.
- X Press the  $\square$  and  $\square$  keys to adjust the absorbance or %transmittance to the desired value.
- X Press the **Set A** soft key or the **Set %T** soft key to set the value and return to the main display.

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### Printing your results

Pressing **Print** sends data to the selected printer. The format of the data is selected using the **Printer Setup** option in the SETUP function. If you press **Print** when the screen is not displaying data, a list of parameters for that screen will be sent to the printer.

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### Using the utility functions

Although you may access the utility functions at any time, you will probably set up these options once and then leave them. To access the utility functions, press **Utility**. You can press  $\square$  or  $\square$  to step through the utilities and select the function you want. Press **Esc** to exit from the Utility menu.

The next few pages list the utility functions that can appear and describe how to use them.

### Concentration units

This utility allows you to select the units of measure used to report concentration values. This utility is available only for concentration mode.

- X Press  $\square$  or  $\square$  until the units screen appears.
- X Press **Change** to select different units.
- X Press  $\square$  or  $\square$  to display the appropriate units.
- X Press **Accept** to select the units or **Esc** twice to exit the Utility menu.

#### Notes:

- X *Available concentration units include: C, ppm, ppb, g/L, mg/L, mg/dL, mg/mL,  $\mu$ g/L, U/L, M, mM,  $\mu$ M and IU*
- X *You can press **PRINT** to obtain a printout of the units available.*

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## Language

This utility allows you to select the language used for screen displays.

- X Press **Change** to view the available languages.
- X Press  or  to step through the list and display the appropriate language.
- X Press **Accept** to select it or **Esc** to exit the Change menu.

## RS232

This utility allows you to set the following RS232 parameters:

- X Baud rate - 300, 600, 1200, 2400, 4800, 9600, 19200
- X Data bits - 8 or 7
- X Parity - off, even, odd
- X Handshake - off, CTS/RTS, XON/XOFF
- X Terminator - CR, LF, CR/LF

- X Press **Change** to select RS232 parameters.
- X Press  or  to display the appropriate values.
- X Press **Accept** to select it or **Esc** to exit the Change menu.
- X Continue for all parameters.

**Note:** You can press **PRINT** to obtain a printout of the RS232 settings.

## Paper Feed

This utility allows you to feed the paper into the printer and to advance the paper in the printer.

- X Press **Feed** to advance the paper in the printer one line at a time.
- X Press **Esc** to exit the Utility menu.

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#### Printer setup

This utility allows you to select the format of the data that will be sent to the printer and to select the printer where the data is sent.


- X Press **Change** to select the printer and printer format.
- X Press  or  to display the printers (OFF, Internal, Parallel, Serial).
- X Press **Accept** to select the printer.
- X Press  or  to display the printer formats (Print WL - Yes/No, Print Units - Yes/No, Decimal Point - . or ,).
- X Press **Accept** to select each printer format or **Esc** to return to the other utilities.

**Note:** You can press **PRINT** to obtain a printout of the printer format settings.

#### Lamp alignment

This utility helps you align the lamp.

- X Press **Align** to display an energy graph like the one below.
- X Adjust the alignment as described in the *Maintenance* section (page 44) until the graph displays a peak energy level.
- X Press **Esc** to exit the lamp alignment utility.

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#### Hours of lamp operation

This utility displays the total hours that the lamp has operated.

- X When you select this utility, a screen appears indicating the number of hours that the lamp in the instrument has been operating (see below).
- X This utility also allows you to reset the counter when you install a new lamp.
  
- X Press **Reset** to set the number of hours to zero.
- X Press **Confirm** and the instrument resets the counter to zero.
- X Press **Esc** to exit the Utility menu or press  or

Lamp hours = 538

Esc      ↓      Reset

to scroll through the other utilities.

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#### Initial wavelength

This utility allows you to set the initial power-on wavelength.

- X When you select this utility, a screen appears indicating the current initial power-on wavelength.
- X Press **Change** to select the power-on wavelength.
- X Press  or  to change the wavelength value.
- X Press **Accept** to accept it, then **Esc** to exit the change menu.

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#### Beeper control

This utility allows you to disable the error beeper.

- X Press **Change** to turn the beeper on or off.
- X Press **Esc** to exit the Utility menu or press  or  to scroll through the other utilities.

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### Standby

This utility allows you to turn off the lamp when there is no activity after a set time. Options for this utility are OFF, 30, 60, 90, 120 150 or 180 minutes.

- X Press **Change** to change the standby time.
- X Press  or  to select the appropriate standby time.
- X Press **Accept** to accept it, then  or  to scroll through the other utilities or **Esc** to exit the Utility menu.

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
### Software revision

This utility allows you to view the software revision of programmed parts in the instrument. When you select **Software revision** from the Utility menu, the software revision is displayed.

- X Press **Esc** to exit the Utility menu or press  or  to scroll through the other utilities.

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
**Accessories**

| <b><u>Part #</u></b> | <b><u>Description</u></b>  |
|----------------------|--|
| 4010                 | Single Cell Holder - holds one 10mm square cuvette, 10mm round test tube or 2 $\equiv$ test tube up to 100mm tall (standard with instrument) |
| 4011                 | Cylindrical Longpath Cell Holder - holds one cylindrical longpath cell up to 50mm long and 22mm in diameter                                  |
| 4012                 | Rectangular Longpath Cell Holder - holds one rectangular longpath cell up to 50mm long and 12mm wide   |
| 4014                 | COD Vial Cell Holder - holds one COD vial up to 104mm tall and 18mm in diameter  |
| 4015                 | 1" Test Tube Holder - holds one 1" pathlength test tube, up to 104mm tall  |
| 4028                 | Filter Holder - holds one filter/lens up to 8mm thick, with dimensions up to 44mm W x 104mm H  |
| 4088                 | 20-Column Internal Printer (user-installable); includes 5 rolls of paper   |
| 333150               | SPECTRONIC Standards   |
| 336041               | Interface cable to connect 20 GENESYS Series to IBM PC/XT computer   |
| 336042               | Interface cable to connect 20 GENESYS Series to IBM/AT computer  |
| 336043               | Interface cable to connect 20 GENESYS Series to non-IBM equipment  |

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

#### Removing the cell holder

- X Firmly grasp the tab at the back of the cell holder and pull forward and up to remove it.

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

#### Installing the cell holder

- X Press down firmly on the cuvette holder until you hear the cell holder snap into place.

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



## Maintenance & Troubleshooting

Your 20 GENESYS<sup>9</sup> spectrometer is durable and reliable, so you can expect trouble-free operation and minimal maintenance. However, periodically, you will need to replace the lamp. You may also need to change the fuse.

This section describes routine maintenance procedures and provides a list of messages generated by the instrument to inform you of errors.

▲ Do not remove the cover of the instrument!

### Routine care

Routine care for your 20 GENESYS<sup>9</sup> spectrometer does not require a lot of time. To help minimize maintenance time and to increase the life and performance of your instrument, please follow these guidelines:

- X Do not use or store the instrument in a corrosive environment.
- X Gently wipe the outside of the instrument with a soft cloth to remove any dust or spills. Water, isopropyl alcohol and other common laboratory cleaning agents may be used if necessary.
- X Always clean up spills immediately to prevent or minimize damage to the instrument.
- X Use water, isopropyl alcohol or other common laboratory cleaning agents to clean the keyboard. It is recommended that you clean spills off the keyboard as soon as they occur.

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
## Changing the lamp

<sup>A</sup>  
The lamp gets very hot during operation. Before removing the lamp, turn off the instrument and allow the lamp to cool down for 10 minutes.

▲ Turn off and unplug the instrument before removing the lamp door.

The instrument and printer contain static-sensitive devices. Take proper electrostatic discharge (ESD) precautions when installing the printer (such as using an ESD wrist strap).

1. **Turn off and unplug the instrument.**
2. Allow the lamp to cool for at least 10 minutes.
3. Use a coin to loosen the door of the lamp compartment, then pull the door up and lift it to remove it.
4. Loosen and remove the thumbscrew (#1, below).
5. Pull the tab (#2) at the top of the lamp holder assembly to lift the assembly out of the instrument.  
**Note:** *The cable remains connected to the instrument.*
6. Pull the lamp straight out to remove it from its socket as shown by the arrow in the figure.

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

▲ Do not touch the new lamp with your fingers!

- 17 Align the two legs protruding from the clean, new lamp with the two holes in the socket and press the lamp securely into the socket (the legs should fit snugly in the holes).
- 18 Use a soft, lint-free cloth to wipe off any fingerprints on the lamp.
- 19 Position the lamp holder assembly pressing the spring (#3) on the right and aligning the pin (#4) with the hole in the assembly.
- 20 Gently slide the assembly into place and replace the thumbscrew, leaving the thumbscrew slightly loose.

▲ The lamp gets very hot during operation.

- 21 Plug the instrument into an appropriate wall outlet and turn on the power.  
**Note:** *If you get an initialization error when you try to start the instrument, turn off and unplug the instrument, then turn the alignment screw  turn in either direction and try again.*
- 22 Press **SETUP** to access the SETUP functions.
- 23 Press  until **Lamp Alignment** appears.
- 24 Press **Align** to display an energy graph like the one below.

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
### Changing the fuse

The fuse is located in the power entry module on the


The instrument fuse must be replaced with the same type and rating fuse. Use IEC standard 5x20mm, Type F fuse, 2.5A, 250VAC, 1500A interrupt rating [Littlefuse #216025].

back of the instrument.

- 86 **Turn off and unplug the instrument.**
- 87 Position the instrument so you can access the power entry module on the back of the instrument.
- 88 Remove the power cord.

If the fuse fails repeatedly, it may indicate a serious problem with the instrument. Contact your service representative as soon as possible.

- 89 Insert a flat-blade screwdriver into the notch on the fuse cover and pry off the cover as shown below.

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



- 103 Unsnap both fuses to remove them (#1, below).
- 104 Insert the new fuses, pushing them in so they snap into place.
- 105 Replace the fuse cover.
- 106 Replace the power cord.
- 107 Plug the instrument back into the appropriate wall outlet and turn on the power.

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
|      |  |
|------|--|
| 4010 | Single cell holder for 10mm square cuvette, 10mm test tube or □□ test tube |
| 4062 | Fuses, package of 4  |
| 4054 | Paper for 4088 Internal Printer, package of 5 rolls                        |

|            |   |
|------------|---|
| 4055       | Long-life paper for 4088 Internal Printer, package of 5 rolls |
| 4001-10001 | 20 GENESYS Spectrometer Operator's Manual                     |
| 4001-10020 | 20 GENESYS Spectrometer Service Manual                        |

#### Messages

This section lists the messages generated to alert you of errors or other abnormal conditions. The instrument recognizes two types of errors. With the first type, the instrument is still functional; with the second, the instrument is not functional until the condition is resolved.

#### *Flashing Data Display*

This condition indicates that the sample has an absorbance or a transmittance value below or above the photometric range of the instrument. The display flashes until the condition is resolved.

#### *Lamp Failure*

This condition indicates that the lamp has failed. The message remains on the display until the lamp is replaced.

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#### *Sample too dark*

This condition indicates that the instrument has been asked to zero a sample with a high absorbance at a low energy point. The instrument beeps three times to announce the message, the message remains on the display for two seconds, then the normal display returns.

#### *Sample too bright*

This condition indicates that the instrument has been asked to zero a sample while the door of the sample compartment is open. The instrument beeps three times to announce the message, the message remains on the display for two seconds, then the normal display returns.

#### *Key not active*

This condition indicates that a non-functional key was pressed. For example, this message would appear if you press **MODE** when viewing the SETUP menu. The instrument beeps three times to announce the message, the message remains on the display for two seconds, then the normal display returns.

#### *Parity error*

This message appears when you try to transmit information from the computer to the instrument and indicates that the setting selected for the parity is incorrect. Use the RS232 utility (see page 26) to select the correct setting.

#### *Framing error*

This message appears when you try to transmit information from the computer to the instrument and indicates that the number of data bits, number of stop bits or both settings are incorrect. Use the RS232 utility to select the correct settings.

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#### *Filter Whl Init Fail/Insufficient Energy*

This message indicates that there is insufficient energy available to initialize the filter wheel during start-up. This can occur when the door of the sample compartment is open during start-up, when a high-absorbing sample is in the sample compartment or when a new lamp has been installed. Turn off and unplug the instrument, turn the alignment screw □ in either direction, then plug in and turn on the instrument again.

- X If the start-up sequence completes successfully, press **SETUP** to access the SETUP functions, then align the lamp.
- X If the start-up sequence does not complete successfully, contact your authorized service representative.

#### *Mono Init Fail/Insufficient Energy*

This message indicates that there is insufficient energy available to initialize the monochromator during start-up. This can occur when the door of the sample compartment is open during start-up or when a new lamp has been installed. Turn off and unplug the instrument, turn the alignment screw □ in either direction, then plug in and turn on the instrument again.

- X If the start-up sequence completes successfully, press **SETUP** to access the SETUP functions, then align the lamp.
- X If the start-up sequence does not complete successfully, contact your authorized service representative.

#### *Hardware failure 13-16*

These conditions occur if the door of the sample compartment is open during the autozero procedure. If one of these messages appears when the door of the sample compartment is closed, contact your authorized service representative.

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